

Material No: 9095.1000
Batch No: 2534605853
Manufactured Date: 2025/12/12
Expiration Date: 2028/12/12
Release Date: 2025/12/15
Revision No.: 0

Certificate of Analysis

Test	Specification	Result
Assay (by GC) (corrected for water)	$\geq 99.7 \%$	100.0
Residue after Evaporation (in ppm)	≤ 2	< 2
Water (H ₂ O)	$\leq 0.03 \%$	0.03
Ultraviolet Absorbance (1.00-cm path vs water):		
at 225 nm	≤ 0.16	0.11
at 254 nm	≤ 0.02	0.01
at 280 nm	≤ 0.01	< 0.01
at 350 nm	≤ 0.01	< 0.01
UV Cut-off, nm	≤ 205	202

For Laboratory, Research or Manufacturing use.
For use in High Performance Liquid Chromatography.
Filtered through a 0.2 micron filter. Packaged under Nitrogen.
Country of Origin: NETHERLANDS
Packaging Site: Gliwice Mfg & DC



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Certificate of Analysis, please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000