

Water
HPLC
For use in Liquid Chromatography (HPLC & UHPLC) &
Spectrophotometry



Material No.: 4218-P04
Batch No.: P24C100302
Manufactured Date: 2024-03-09
Expiration Date: 2028-03-08
Revision No.: 0

Certificate of Analysis

Test	Specification	Result
UV Absorbance on Gradient Elution – Larg	≤ 0.001	≤ 0.001
UV Absorbance on Gradient Elution – Larg		0.003
Fluorescence Trace Impurities, measured	≤ 0.1 ppb	≤ 0.1 ppb
Fluorescence Trace Impurities, measured	≤ 0.2 ppb	≤ 0.2 ppb
Filtration Test	Passes test	Passes test
pH, For Information Only		6.4
Residue after Evaporation	≤ 2.0 ppm	0.8 ppm

For Laboratory, Research, or Manufacturing Use

Country of Origin: India

Girish Prajapati
Site Quality Manager

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA 19087. U.S.A. Phone 610.386.1700

Page 1 of 1