

Petroleum Ether, 35–60°C.  
BAKER ANALYZED® A.C.S. Reagent



Material No.: 9268-05  
Batch No.: 24K2561019  
Manufactured Date: 2024-11-22  
Retest Date: 2029-11-21  
Revision No.: 0

## Certificate of Analysis

Meets ACS Reagent Chemical Requirements

| Test                            | Specification  | Result      |
|---------------------------------|----------------|-------------|
| ACS – Color (APHA)              | $\leq 10$      | <5          |
| ACS – Residue after Evaporation | $\leq 0.001\%$ | <0.001 %    |
| ACS – Acidity                   | Passes Test    | Passes Test |
| ACS – Boiling Range(Lower)      | 35 – 60 °C     | 38 °C       |
| ACS – Boiling Range(Upper)      | 35 – 60 °C     | 59 °C       |

For Laboratory, Research, or Manufacturing Use

Country of Origin: United States  
Packaging Site: Paris Mfg Ctr & DC

  
Michelle Bales  
Sr. Manager, Quality Assurance

For questions on this Certificate of Analysis please contact Technical Services at 855.282.6867 or +1.610.386.1700

Avantor Performance Materials, LLC

100 Matsonford Rd, Suite 200, Radnor, PA, 19087, U.S.A. Phone 610.386.1700